0 Sandersweg 15 in 21680 Stade 2 04141- 535 9810 ^ 04141 - 5359830 S www.johannisheim-stade.de

Johannisheim e.V. - Sandersweg 15 in 21680 Stade

BEO MedConsulting Berlin GmbH attn: Mr. Friedrich Helmholtzstr. 2 Nursing services Sandra Idem 2 04141-5359828 El sandra.idem@johannisheim-stade.de

05/12/2018

10587 Berlin - Charlottenburg

Evidence of the medical benefit of the "Ride Forward Cushion" anti-decubitus foam seat cushion to assist in the prevention and treatment of pressure ulcers by Ride Designs, USA

Dear Mr. Friedrich,

I am pleased to send you today my statement on the proof of the medical benefit of the "Ride Forward Cushion" anti-decubitus seat cushion. The test specimen was CE-marked and provided with instructions for use. After the test period of several weeks, the nursing college and I arrived at the results presented in the following.

Should you have any questions about the results presented, please contact me.

I hope I can help you with my explanations.

Best regards,

Sandra Idem

- Nurse & nursing services -

1. The "Ride Forward Cushion" anti-decubitus seat cushion incl. protective cover was evaluated. The product is intended to assist in the prevention and treatment of pressure ulcers in accordance with the manufacturer's product documentation. It is offered in the following dimensions:

Ride Forward Cushion (25 x 25 cm) Art. no.: FCD-1010 Ride Forward Cushion (25 x 30 cm) Art. no.: FCD-1012 Ride Forward Cushion (30 x 30 cm) Art. no.: FCD-1212 Ride Forward Cushion (30 x 36 cm) Art. no.: FCD-1214 Ride Forward Cushion (30 x 41 cm) Art. no.: FCD-1216 Ride Forward Cushion (36 x 36 cm) Art. no.: FCD-1414 Ride Forward Cushion (36 x 41 cm) Art. no.: FCD-1416 Ride Forward Cushion (36 x 46 cm) Art. no.: FCD-1418 Ride Forward Cushion (38 x 38 cm) Art. no.: FCD-1515 Ride Forward Cushion (38 x 43 cm) Art. no.: FCD-1517 Ride Forward Cushion (41 x 41 cm) Art. no.: FCD-1616 Ride Forward Cushion (41 x 46 cm) Art. no.: FCD-1618 Ride Forward Cushion (41 x 51 cm) Art. no.: FCD-1620 Ride Forward Cushion (43 x 43 cm) Art. no.: FCD-1717 Ride Forward Cushion (46 x 41 cm) Art. no.: FCD-1816 Ride Forward Cushion (46 x 46 cm) Art. no.: FCD-1818 Ride Forward Cushion (46 x 51 cm) Art. no.: FCD-1820 Ride Forward Cushion (51 x 41 cm) Art. no.: FCD-2016 Ride Forward Cushion (51 x 46 cm) Art. no.: FCD-2018 Ride Forward Cushion (51 x 51 cm) Art. no.: FCD-2020 incl. inner and outer protective cover

The scope of delivery includes the cushion with an inner and outer protective cover.

2. Information on the tester and institution

The Ride Forward Cushion anti-decubitus seat cushion was tested at the Johannisheim care home in Stade. The facility provides individual inpatient care for the elderly and the disabled (at all levels of care) in various living areas. The individual wishes and needs of guests and residents are at the heart of the nursing care.

Ms. Sandra Idem acts as the expert. She is a registered nurse and provider of nursing services. Ms Idem has been testing anti-decubitus systems for many years.

3. Information on the period of use/observation

The anti-decubitus seat cushion was tested from 15 August 2018 to 8 November 2018 in our facility with six residents over a period of six weeks each.

4. Description of test subjects and application history

Subjects were between 65 and 91 years of age and had different underlying medical conditions. These were predominantly typical, predominantly age-related neurological, internal and neurodegenerative diseases. The risk of pressure ulcers was assessed according to the national expert standard for pressure ulcer prevention based on nursing expertise. In addition, the Braden scale was used as a standardised assessment.

Residents	01	02	03
Test location	Stade		
Test period	from 15 August 2018 to 26 Septemb	per 2018	
Gender	Male	Female	Female
Age	76	65	67
Height in cm	170	178	170
Weight in kg	98	82	107
Diagnoses, specia features of the card situation	 eregion Immobility syndrome Transtibial amputation left Forefoot amputation right Diabetes mellitus type II 		 ACID occlusion and lysis therapy with right middle cerebral artery infarction right. Alcohol abuse Art. hypertension Extensive middle cerebral artery infarction with paresis left Slipped disc 3/2011 Depression Diab. mellitus type II Facial nerve paresis left Scheuermann's disease Nicotine abuse
Level of care	• IV	• V	• IV

Pressure ulcer risk assessment through nursing expertise	High risk	High risk	High risk
Braden-scale assessment of pressure ulcer risk	Medium risk (14 points)	Very high risk (07 points)	Medium risk (14 points)
Special risk factors	 Diabetes mellitus Bed confinement with very severe movement restrictions Paralysis Limited sensory faculties No conscious relief of the sitting area Pelvic obliquity due to amputation 	 Danger of malnutrition Bed confinement with very severe movement restrictions Presence of contractures Limited sensory faculties No conscious relief of the sitting area There is no stability while sitting Limited balance Inadequate change of sitting position 	 Diabetes mellitus Very pronounced movement restrictions Paralysis and spasticity Limited sensory faculties Pain symptoms No conscious relief of the sitting area The resident can change the sitting position at times
Length of use of test specimen per day	 Continuously for 4 hours at a time, 8 hours daily in total 	Continuously for 2 hours,2 hours daily in total	Continuously for 10 hours,10 hours daily in total
Pressure ulcer present	Yes	Yes	Yes
Category according to EPUAP		I buttock region	II buttock region
Healing progress	Positive	Positive	Positive
Has there been a pressure ulcer in the past?	No	No	No
Estimation of pressure relief with the aid	Good - no redness or maceration has occurred since use of the cushion	Good - no unwanted skin changes have occurred	Good - no unwanted skin changes have occurred
Positioning	• The resident can be positioned in the usual way	• The resident can be positioned in the usual way	• The resident can be positioned in the usual way
Important remarks	• None	• None	• None

Residents	04	05	06
Test location	Stade		
Test period	from 15 August 2018 to 26 September 2018	from 27 September 2018 to 8 November 2018	from 27 September 2018 to 8 November 2018
Gender	Male	Male	Male
Age	91	73	80
Height in cm	175	168	198
Weight in kg	79	86	109
Diagnoses, special	 Post apoplexy with left 	Art. hypertension	•
features of the care situation		 Art. hypertension Diab. Mellitus Seizures Left middle cerebral artery infarction with right hemiparesis 	
Level of care	• IV	• IV	• IV
	High risk	Medium risk	Medium risk
Braden-scale assessment of pressure ulcer risk	General risk (15 points)	Low (18 points)	Medium risk (14 points)
Special risk factors	 Resident is advanced in years Moderate movement restrictions Paralysis Limited sensory faculties Partial conscious relief of the sitting area 	 Moderate movement restrictions Paralysis and contractures Limited sensory faculties 	 Resident is advanced in years Diabetes mellitus Limited sensory faculties No conscious relief of the sitting area Asymmetrical sitting
Length of use of test specimen per day	 Continuously for 4 hours at a time, 4 hours daily in total 	• 8 hours daily in total	 Continuously for 4 hours, 8 hours daily in total
Pressure ulcer present		Yes	No
Category according to	Category I, sacral region	Category I, sacral region	Not applicable
EPUAP Healing progress	Positive	Positive	Not applicable
Has there been a pressure ulcer in the past?	No	No	No
Estimation of pressure relief with the aid	Good - no redness or maceration has occurred since use of the cushion	Good - no unwanted skin changes have occurred	Good - no unwanted skin changes have occurred
Positioning	• The resident can be positioned in the usual way	Resident can move into position independently	• The resident can be positioned in the usual way

5. Information on the test equipment

The Ride Forward Cushion anti-decubitus seat cushion was tested at our facility in Stade.

Our facility

Our care team at Johannisheim nurtures and cares for people with mental and physical complaints at all levels of care. Johannisheim has capacity for 124 residents. 90 single and 17 double rooms of different room size(s) are available. Our services cover short and long term care and inclusive dementia care. Johannisheim also cares for residents who have suffered from a stroke, Parkinson's disease, apallic conditions, Alzheimer's disease and MS, etc.

Our philosophy

Johannisheim is a social welfare institution based on the Christian ideology. This affects the basic social diaconal attitude and actions in all areas of the facility's life and work.

The job is a vocation for us

Everyone has an unrestricted right to a dignified life - even at an advanced age. This conviction determines how we deliver our service to people in care and support. We want to contribute to humanity in our society through expert support and personal care - by helping and listening.

Johannisheim sees itself as a "company of humanity"

- Our work is characterized by appreciation for people in need of care.
- We want to combine our professional competence with personal care.
- We see our cooperation as an important asset for our residents.
- The satisfaction of the people entrusted to us and who we work to help is an important benchmark for us.
- Good cooperation is an essential foundation for our welfare services.

Feedback is a development aid for us

The senior advisory board freely chosen in our company is also our in-house supervisory body and monitors all services and activities in the business, such as care, catering, supervision, accommodation, care contracts, construction measures, investments, projects, etc., but also the socially responsible treatment of the residents entrusted to us. Our senior advisory board, as well as relatives, guests and visitors, improve our quality of work and enrich our development through critical feedback of a positive and negative nature.

6. Statement on the effectiveness of the test specimen

The test specimen was used in our facility to support the management of pressure ulcers for six residents.

The risk of pressure ulcers is assessed in our facilities based on professional nursing expertise. In line with this assessment of the risk of pressure ulcers, 4 residents posed a "high risk" and 2 a "medium risk".

As part of this test, the Braden scale was also used to determine the risk of pressure ulcers. In accordance with the Braden scale, all residents in our facility were at risk of developing a pressure ulcer. The following risks were documented in accordance with Braden:

- 1 resident had a "general risk"
- 1 resident had a "low risk"
- 3 residents had a "medium risk"
- 1 resident had a "very high risk"

5 residents also had a pressure ulcer in the sacral region - 3 category I (EPUAP) and 2 category II (EPUAP) pressure ulcers.

In addition to using the test specimen, the residents were positioned at the necessary time intervals and received individual wound care, as well as specific skin care.

Subject	Pressure ulcer localisation	Wound initial situation	Wound situation after the three-week test period
1	Posterior region	Category II EPUAP, fibrinous,	Category II EPUAP,
	region	no redness (depth x length x	granulating, no redness (depth
		width in cm) 0.2 x 2.0 x 1.5	x length x width in cm) 0.1 x
		Category I EPUAP,	Category I EPUAP,
2	Posterior region	(depth x length x width in cm)	(depth x length x width in cm)
		0.0 x 1.0 x 2.0	0.0 x 1.0 x 1.0
		Category II EPUAP, fibrinous,	Category II EPUAP, fibrinous,
		signs of wound infection (depth	signs of wound infection
3	Posterior region	x length x width in cm) 0.1 x	(depth x length x width in cm)
		2.0 x 1.0	0.1 x 1.6 x 0.8
		Category I EPUAP, (depth x	Healed
4	Posterior region	length x width in cm) 0.0 x 2.0	
		x 2.0	
		Category I EPUAP, (depth x	Healed
5	Posterior region	length x width in cm) 0.0 x	
		2.0 x 2.0	

The table presentation of the wound progression shows that a visible improvement in the wound condition was observed in the case of pressure ulcers or two of the pressure ulcers healed.

Important components of our facility's pressure ulcer management include not only local wound treatment, but also all preventative measures such as nutrition, skin care, promoting movement, pressure distribution through individual positioning and the use of anti-decubitus mattresses and seat cushions. Consequently, the use of an anti-decubitus seat cushion is an important part of the overall nursing care concept.

From the observations and results of the test, it must be summarised that the nursing team is satisfied with the medical benefit of the test specimen, as no pressure ulcers or recurrences have occurred and the treatment of pressure ulcers has been positively supported.

7. Assessment of the benefits against the risks

At our facility, we were not able to observe any hazards arising from the product, either for patients or care staff.

The manufacturer's instructions for use are easily understandable and contain the most important information about product properties and materials used, as well as instructions for cleaning, maintenance and care.

8. Evaluation of desired and undesirable consequences

There were no undesirable consequences when using the product. Furthermore, for example, no negative changes in the microclimate were observed due to the occurrence of specific moisture on the skin, although, during the test (in August and September), there were sometimes high summer outside temperatures.

9. **Product description**

The "Ride Forward Cushion" anti-decubitus wheelchair seat cushion is a medical device intended for the prevention, treatment or alleviation of pressure ulcers in the sacral region.

The cushion consists of:

- The actual seat cushion core
- An inner cover
- An outer cover.

The seat cushion core consists of two layers. The upper part is made of viscoelastic foam, the lower part of polyurethane foam. The foam core is anatomically preformed.

The seat cushion core is covered by a first cover made of polyether. It is equipped with a zipper and protects the core from the penetration of liquids.

The outer cover, which is also removable, is made of air-permeable fabrics that are intended to counteract the increased occurrence of moisture on the skin. The underside is also slip-resistant, preventing the cushion from "slipping" on the support (e.g. the seat of a wheelchair).

9.1 Changing covers

The construction of the two covers allows them to be changed as they both have a zipper.

9.2 Adaptability of the product

The product cannot be adjusted by the user.

To use the cushion, simply place it on the seat base. The "Ride Forward Cushion" anti-decubitus seat cushion has a maximum user weight. It can be used for patients weighing 5 to 113 kg.

9.3 Usability of adjustable values by the user The product cannot be

changed by the user.

9.4 Microclimate characteristics

There was no increased occurrence of moisture in the area of the seat surface despite high temperatures.

9.5 Cleaning and disinfection

It is possible to clean and disinfect the cushion as well as the covers (see instructions for use).

9.6 Assessment of the pressure relief

The pressure-distributing effect is evaluated as positive, as no new skin irritation was observed in any of the at-risk patients during the test period for the test specimen. The existing pressure ulcers also showed a healing tendency or healed.

9.7 Ability to change positions as required

The usual seated position changes could be carried out.

9.8 Changes in the position support intervals (positioning intervals)

The position support intervals remained unchanged for all subjects during the test phase.

9.9 Estimation of the applied shear forces

The cushion does not adversely affect shear forces.

9.10 Changes in sitting stability

The sitting stability of the test subjects has not changed negatively.

9.11 Changes in sitting symmetry

The sitting symmetry of the test subjects has not changed negatively.

10. Overall assessment

The "Ride Forward Cushion" anti-decubitus seat cushion was successfully tested in our care facilities. Our residents were satisfied with the product and the resulting sitting quality. The nursing specialists are positive towards the product after several weeks of testing. Handling is easy to assess for both users and carers.

The test specimen was used by high-risk patients over a period of six weeks. Five subjects had a category I and II (EPUAP) pressure ulcer in the sacral region. Progress in wound healing was observed in all of them. Two of the pressure ulcers healed.

Throughout the period of use, none of our residents experienced new pressure-related skin lesions, which proved the effectiveness of the cushion.

Based on the results and experience gained during use of the cushion being tested, we assess the product as suitable for use in pressure ulcer management (pressure ulcer prevention and treatment) in accordance with the manufacturer's indications. The use of the cushion can be advocated for both domestic and inpatient use.

The content of this report or parts thereof may not be published, made accessible to the public or a group of people outside the company or used for advertising purposes without the written consent of the author.

Stade, 04/12/2018

Verification of medical benefit / for products in product group 11

Page 1 of 1

beo berlin

Test plan for producing a report on the medical benefit

for

Product designation: Anti-decubitus foam seat cushion

Type/version:	Ride Forward Cushion (25 x 25 cm)	Art. no.:	FCD-1010
	Ride Forward Cushion (25 x 30 cm)	Art. no.:	FCD-1012
	Ride Forward Cushion (30 x 30 cm)	Art. no.:	FCD-1212
	Ride Forward Cushion (30 x 36 cm)	Art. no.:	FCD-1214
	Ride Forward Cushion (30 x 41 cm)	Art. no.:	FCD-1216
	Ride Forward Cushion (36 x 36 cm)	Art. no.:	FCD-1414
	Ride Forward Cushion (36 x 41 cm)	Art. no.:	FCD-1416
	Ride Forward Cushion (36 x 46 cm)	Art. no.:	FCD-1418
	Ride Forward Cushion (38 x 38 cm)	Art. no.:	FCD-1515
	Ride Forward Cushion (38 x 43 cm)	Art. no.:	FCD-1517
	Ride Forward Cushion (41 x 41 cm)	Art. no.:	FCD-1616
	Ride Forward Cushion (41 x 46 cm)	Art. no.:	FCD-1618
	Ride Forward Cushion (41 x 51 cm)	Art. no.:	FCD-1620
	Ride Forward Cushion (43 x 43 cm)	Art. no.:	FCD-1717
	Ride Forward Cushion (46 x 41 cm)	Art. no.:	FCD-1816
	Ride Forward Cushion (46 x 46 cm)	Art. no.:	FCD-1818
	Ride Forward Cushion (46 x 51 cm)	Art. no.:	FCD-1820
	Ride Forward Cushion (51 x 41 cm)	Art. no.:	FCD-2016
	Ride Forward Cushion (51 x 46 cm)	Art. no.:	FCD-2018
	Ride Forward Cushion (51 x 51 cm)	Art. no.:	FCD-2020

<u>incl. cover</u>

Tested accessories: None

Manufacturer:	Ride Designs a branch of Aspen Seating, LLC, 8100 SouthPark Way, Bldg. C-400, Littleton, Colorado 80120 USA
Client:	BEO MedConsulting Berlin GmbH, Helmholtzstr. 2, 10587 Berlin

Reviewed by: Johannisheim, Sandersweg 15, 21680 Stade

Please create a written report with the following content:

- 1. Information on product designation, type/design with item no., the tested accessories where applicable with item no., manufacturer and customer, e.g. as listed above.
- 2. Brief description of your facility/institution that describes the independence and reputation of your institution and information on the tester.
- 3. Information on the period of use/observation *Note: A minimum observation period of 3 weeks per patient is recommended.*
- 4. Details of the users/subjects used for the test and the number of subjects. *Note: The test must be carried out taking into account the intended user group. We recommend including at least 5 patients in the observations.*
- 5. Information on the test environment (e.g. clinic, home setting) Note: The studies/tests must be carried out under general living conditions/home environment or it must be possible to apply the results to these areas.
- 6. Statement on the entire indication area specified by the manufacturer and the effectiveness of the product in combination with all tested components (e.g. covers, accessories, etc.), with:
 - 6.1. Description of the intended goal and
 - Ensuring medical treatment (decubitus therapy) and/or
 - Prevention of a disability (decubitus prophylaxis) and/or
 - Treatment of a disability (e.g. improving freedom of movement by changing position)
 - 6.2. Clear description of the goal(s) achieved and outlined under 6.1., along with reason why the tester believes that the product is suitable for use **in the home environment** with the indications made by the manufacturer, complete with clinical endpoints, the pressure ulcer incidences (occurrence of pressure ulcer stage II or higher) to verify the effectiveness at preventing pressure ulcers, or, for products intended for pressure ulcer treatment, the healing rate (specifying the measurement method use) of an occurring pressure ulcer. Here, the indication area used must be clearly formulated and outlined based on standardised risk scales the Braden scale in the case of geriatric patients and other approved scales for other patient groups; if the product is intended for pressure ulcer treatment, a description of the relevant pressure ulcer stages must also be provided in accordance with EPUAP.
- 7. An assessment of the benefits against the risks.
- 8. An assessment of the desired and undesirable consequences ("outcomes")
- 9. Description of the product with statement on the following points:
 - 9.1. Can the cover be changed by the user?
 - 9.2. Can the product be adapted to the individual weight of the patient and the relevant load situation (manually or automatically)? Or can the product be selected based on the weight of the patient?
 - 9.3. Can all individually adjustable values be operated by the user?
 - 9.4. Observable microclimatic properties
 - 9.5. Further information on the handling of the product where applicable, e.g. cleaning/disinfection, comfort for the user/patient and care staff.
 - 9.6. Assessment of the pressure relief
 - 9.7. Ability to change positions as required
 - 9.8. Changing the position support intervals (positioning intervals)
 - 9.9. Estimation of the applied shear forces
 - 9.10. Changes in sitting stability
 - 9.11. Changes in sitting symmetry

10. Overall assessment

The report can be written in German or English and pictures and evaluation tables may be added where applicable.

Confirmation of the test facility /expert

I hereby confirm that the assessment of the medical benefit for the product listed above was created according to the test plan provided.

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