Centre for ventilation and intensive care in "Storkower Bogen" GmbH Franz-Jacob-Straße 4D ♦ 10369 Berlin

BEO MedConsulting Berlin GmbH attn: Mr. Friedrich Helmholtzstr. 2

10587 Berlin - Charlottenburg

Sebastian Kruschwitz Specialist Coordinator Wound Management, Deputy Case Manager

Tel 030 857 4725 34 Fax 030 857 4725 19

E-mail: sebastian.kruschwitz@zentrum-berlin.com

www.zbi-gruppe.com

Thursday, 6 June 2019

Evidence of the medical benefit of the "Ride Java Cushion" anti-decubitus foam seat cushion to assist in the prevention and treatment of pressure ulcers by Ride Designs, USA

Dear Mr. Friedrich,

I am pleased to send you today my statement on the proof of the medical benefit of the "Ride Java Cushion" anti-decubitus seat cushion. The test specimen was CE-marked and provided with instructions for use. After the test period of several months, the nursing college and I arrived at the results presented in the following.

Should you have any questions about the results presented, please contact me. I hope I can help you with my explanations.

Best regards,

Sebastian Kruschwitz

- Head of Wound Management Department - Kind regards,

1. The "Ride Java Cushion" anti-decubitus seat cushion incl. protective cover was evaluated. The product is intended to assist in the prevention and treatment of pressure ulcers in accordance with the manufacturer's product documentation. It is offered in the following dimensions:

. no.:	<u>JC-1414</u>
. no.:	JC-1416
. no.:	<u>JC-1515</u>
. no.:	<u>JC-1517</u>
. no.:	<u>JC-1616</u>
. no.:	<u>JC-1618</u>
. no.:	<u>JC-1620</u>
. no.:	<u>JC-1717</u>
. no.:	<u>JC-1816</u>
. no.:	<u>JC-1818</u>
. no.:	<u>JC-1820</u>
. no.:	<u>JC-2016</u>
. no.:	<u>JC-2018</u>
. no.:	<u>JC-2020</u>
_	. no.:

The cushion is offered with a cover.

2. Information on the tester and institution

The "Ride Java Cushion" anti-decubitus seat cushion has been tested at the "Storkower Bogen" ventilation and critical care centre in Berlin.

Since its foundation in 2006, the centre has established itself as a full inpatient and outpatient nursing facility in Berlin. As a solitary facility, it stands for an innovative, rehabilitative and specialist nursing concept for people with outpatient ventilation and intensive care needs, for people in the so-called unresponsive wakefulness syndrome (UWS) and for people in a minimally conscious state (MCS).

The head of the Wound Management Department, Sebastian Kruschwitz, acts as the expert. He is a nurse and wound expert/ICW care specialist.

3. Information on the period of use/observation

The anti-decubitus seat cushion was tested in our facility by 14 patients/residents over a period of six weeks each from 5 November 2018 to 14 February 2019.

4. Description of test subjects and application history

Subjects were between 48 and 83 years of age and had different underlying medical conditions. These were neurological and internal illnesses. Patients were completely limited or bedridden in their mobility.

The risk of pressure ulcers was assessed on the basis of the Braden scale and nursing expertise according to the national expert standard for pressure ulcer prevention. Subjects had different assessments regarding the risk of pressure ulcers. All were at risk of developing a pressure ulcer.

Residents	01	02	03	
Test location		Berlin		
Test period	fro	m 5 November 2018 to 15 February 2	2019	
Gender	Male	Male	Female	
Age	51	62	57	
Height in cm	167	180	168	
Weight in kg	76.6	87	62	
Diagnoses, special features of the care situation	 Mucopolysaccharidosis IV Tetraparesis and plegia Damage to the spinal cord of C6 -C8 with multiple bone irritations Post decubitus III and flap surgery Ventilation Tracheostoma 	 Multiple stroke syndrome Global aphasia Dysphagia Complete paraplegia sub Th4 with spinal ischaemia 	Basal ganglia bleeding, right, with ventricle collapse with arteriovenous malformation Diabetes mellitus type II	
Level of care	• V	• V	• V	
Pressure ulcer risk assessment through nursing expertise	Very high risk	Very high risk	Very high risk	
Braden-scale assessment of pressure ulcer risk	Very high risk (6 points)	Very high risk (6 points)	Very high risk (6 points)	

Special risk factors	There is a risk of malnutrition Permanently low blood pressure Unstable circulation Subject lies on a pressure ulcer healed to a significant degree Bedridden Subject is completely immobile Contractures, spasticity and paralysis are present Sensory perception is limited There is no stability while sitting Limited balance Inadequate change of sitting position Asymmetrical pelvic obliquity Hip joint stiffness Impairment of abduction and adduction	Danger of malnutrition Subject lies on a pressure ulcer healed to a significant degree Bedridden Subject is completely immobile Paralysis and contractures are present Sensory perception is limited There is no stability while sitting Limited balance Inadequate change of sitting position Asymmetrical pelvic obliquity Hip joint stiffness Impairment of abduction and adduction	 Danger of malnutrition Diabetes mellitus Bedridden Subject is completely immobile Sensory perception is limited Subject sits on a pressure ulcer healed to a significant degree There is no stability while sitting Limited balance Inadequate change of sitting position Asymmetrical pelvic obliquity Impairment of abduction and adduction
Length of use of test specimen per day	 Continuously for 2 hours at a time, 4 hours daily in total 	 Continuously for 2 hours, 2 hours daily in total 	 3 to 4 hours continuously, 3 to 4 hours daily in total
Pressure ulcer present	No	Yes	Yes
Category according to EPUAP	Not applicable	III sacral region	II buttock region
Healing progress	Not applicable	Positive	Positive
Has there been a pressure ulcer in the past?	Sacral region category III with plastic cover	Os sacrum category III	No
Estimation of pressure relief with the aid	Achieved	Achieved	Achieved
Positioning	The resident can be positioned in the usual way	The resident can be positioned in the usual way	The resident can be positioned in the usual way
Important remarks	• None	 Patient had a much safer sitting position The positioning interval could be extended from 2 to 3 hours No shear forces in the area of the sacrum as free positioning in this area 	Complete pressure relief around the sacrum The positioning interval could be extended by 1-2 hours

Residents	04	05	06
Test location		Berlin	
Test period	fr	rom 5 November 2018 to 15 February 20	019
Gender	Male	Male	Female
Age	75	76	56
Height in cm	182	190	162
Weight in kg	97	80	69
Diagnoses, special features of the care situation	 Respiratory insufficiency with resp. acidosis Congenital spastic quadriplegia 	 Respiratory insufficiency Chron. renal failure Diabetes mellitus type II 	Respiratory insufficiency COPD Gold -IV
Level of care	• V	• V	• III
Pressure ulcer risk assessment through nursing expertise	Very high risk	Very high risk	Medium risk
Braden-scale assessment of pressure ulcer risk	Very high risk (6 points)	Very high risk (6 points)	General risk (17 points)
Special risk factors	Subject lies on a pressure ulcer healed to a significant degree Bedridden Subject is completely immobile Contractures, spasticity and paralysis are present Sensory perception is limited There is no stability while sitting Limited balance Inadequate change of sitting position Asymmetrical pelvic obliquity Hip joint stiffness Impairment of abduction and adduction	Danger of malnutrition Permanently unstable circulation Diabetes mellitus May lie or sit on a healed, highergrade pressure ulcer Bedridden Subject is completely immobile Sensory perception is limited Inadequate change of sitting position Limited balance	Presence of malnutrition and the risk of malnutrition Unstable circulation Moderate movement restrictions Pain symptoms
Length of use of test specimen per day	 Continuously for 2 hours at a time, 3-4 hours daily in total 	 Continuously for 4-5 hours, 4-5 hours daily in total 	Continuously for 2-3 hours,6 hours daily in total
Pressure ulcer present		No	No
Category according to EPUAP		Not applicable	Not applicable
Healing progress	Not applicable	Not applicable	Not applicable

Has there been a pressure ulcer in the past?		Yes, Os sacrum category III	No
Estimation of pressure relief with the aid	Achieved	Achieved	Achieved
Positioning	*	The resident can be positioned in the usual way	The resident can be positioned in the usual way
Important remarks	_ = =	Positioning intervals The positioning intervals could be extended by 2-3 hours	Resident was able to sit longer Resident is more stable

Residents	07	08	09
Test location	Berlin		<u> </u>
Test period	from 5 November 2018 to 15 February 2019		
Gender	Female	Male	Female
Age	56	65	57
Height in cm	172	172	163
Weight in kg	86	80	81
Diagnoses, special features of the care situation	Trepanatius Defelt le. temporal Decompressive craniectomy due to subdural haematoma	 Hypoxic brain damage after resuscitation Bleeding anaemia requiring transfusion Global respiratory insufficiency COPD Diabetes type 11 	Nosocomial pneumonia, both sides Subarachnoid haemorrhage in M2 aneurysm, right
Level of care	• V	• V	• V
Pressure ulcer risk assessment through nursing expertise	Very high risk	Very high risk	Very high risk
Braden-scale assessment of pressure ulcer risk	Very high risk (6 points)	Very high risk (6 points)	Very high risk (6 points)
Special risk factors	Bedridden Subject is completely immobile Contractures are present There is no stability while sitting Limited balance Inadequate change of sitting position Impairment of abduction and adduction	Diabetes mellitus May lie or sit on a healed, highergrade pressure ulcer Bedridden Subject is completely immobile Paralysis and contractures are present Sensory perception is limited No conscious relief of the sitting area No stability while sitting Inadequate change of sitting position Asymmetrical sitting	May lie or sit on a healed, higher-grade pressure ulcer Bedridden Subject is completely immobile Paralysis present Sensory perception is limited. There is no stability while sitting. Asymmetrical sitting Abduction and adduction are impaired Balance is limited

Length of use of test specimen per day Pressure ulcer	 Continuously for 6 hours at a time, 3-4 hours daily in total 	 Continuously for 2-3 hours, 3 hours per day in total 	 Continuously for 2-3 hours, 3-4 hours daily in total
Category according	Not applicable	IV sacral region	Not applicable
Healing progress Has there been a pressure ulcer in the past?	Positive No	Positive No	Not applicable
Estimation of	Good	Good	Good
Positioning	The resident can be positioned in the usual way	The resident can be positioned in the usual way	The resident can be positioned in the usual way
Important remarks	The positioning interval could be extended	• None	Positioning intervals could be extended

Residents	10	11	12
Test location	Berlin		
Test period	fr	om 5 November 2018 to 15 February 20)19
Gender	Male	Female	Male
Age	48	83	72
Height in cm	185	164	156
Weight in kg	88	79	76
Diagnoses, special features of the care situation	Anoxic brain damage Acute respiratory insufficiency	 Tetraplegia in absolute spinal canal stenosis with myelon compression C 2/3 - C 4/5 Diabetes mellitus type 11 	Hypoxic brain damage Condition after resuscitation
Level of care	• V	• V	• V
Pressure ulcer risk assessment through nursing expertise	Very high risk	Very high risk	Very high risk
Braden-scale assessment of pressure ulcer risk	Very high risk (6 points)	Very high risk (6 points)	Very high risk (6 points)
Special risk factors	 Bedridden with total immobility Contractures are present Sensory perception is limited 	 Resident is advanced in years Diabetes mellitus Danger of malnutrition 	Subject lies or sits on a healed, higher-grade pressure ulcer Bedridden with total immobility

	Limited sensory faculties No conscious relief of the sitting area There is no stability while sitting Balance is limited Inadequate change of sitting position Asymmetrical pelvic level Abduction and adduction are impaired	Subject lies or sits on a healed, higher-grade pressure ulcer Bedridden with total immobility Paralysis and contractures Sensory perception is limited No stability while sitting Limited balance Inadequate change of sitting position Pelvic obliquity Hip joint stiffness Impairment of adduction and abduction	Presence of contractures Pain symptoms Limited sensory perception Limited balance Inadequate change of sitting position Hip joint stiffness Impairment of adduction and abduction
Length of use of test specimen per day	 Continuously for 3 hours at a time, 3-4 hours daily in total 	 Continuously for 2 hours, 2-3 hours daily in total 	Continuously for 2 hours, 2-3 hours daily in total
Pressure ulcer present	No	No	No
Category according to EPUAP	Not applicable	Not applicable	Not applicable
Healing progress	Not applicable	Not applicable	Not applicable
Has there been a pressure ulcer in the past?	No	Yes, formerly Category III sacral region	Yes, formerly Category III sacral region
Estimation of pressure relief with the aid	Good	Good	Good
Positioning	The resident can be positioned in the usual way	The resident can be positioned in the usual way	The resident can be positioned in the usual way
Important remarks	The positioning intervals could be extended	 Extension of mobilisation times Stable sitting No sliding in the wheelchair	Stable sitting - no sliding in the wheelchair Mobilization times could be extended by up to 2 hours

Residents	13	14
Test location		Berlin
Test period	fro	n 5 November 2018 to 15 February 2019
Gender	Female	Male
Age	82	76
Height in cm	162	186
Weight in kg	69	90
	ICB ri. frontoparietal	Hypoxic brain damage
Diagnoses, special	Coma grade III	 Condition after resuscitation
features of the care	Tetraplegia	Respiratory insufficiency
situation	Dysphagia	

	Gastrostomy	
	Tracheostoma	
Level of care	• V	• V
Pressure ulcer risk assessment through nursing expertise	Very high risk	Very high risk
Braden-scale assessment of pressure ulcer risk	Very high risk (6 points)	Very high risk (6 points)
Special risk factors	 Resident is advanced in years There is malnutrition Bedridden with total immobility Signs of paralysis are present Sensory perception is limited Sensory stimuli in the buttocks area cannot be perceived No conscious relief of the sitting area Balance is limited No active change of sitting position 	 Permanently unstable circulation Diabetes mellitus Resident lies on a healed, higher-grade pressure ulcer Bedridden with total immobility Paralysis, contractures and spasticity are present Pain symptoms Sensory perception is limited There is no stability while sitting Balance is limited Pelvic obliquity Stiffening of joints in hips Impairment of adduction and abduction
Length of use of test specimen per day	Continuously for 2-3 hours at a time,4-5 hours daily in total	Continuously for 8 hours,8 hours daily in total
Pressure ulcer present	No	No
Category according to EPUAP	Not applicable	Not applicable
Healing progress	Not applicable	Not applicable
Has there been a pressure ulcer in the past?	No	Category III Os sacrum
Estimation of pressure relief with the aid	Good	Good
Positioning	The resident can be positioned in the usual way	The resident can be positioned in the usual way
Important remarks	Extension of mobilisation phases	 Stable sitting No patient slippage Positioning intervals could be extended

5. Information on the test equipment

The "Ride Java Cushion" anti-decubitus seat cushion was tested at our facility in Berlin.

For specialised nursing care of people with outpatient ventilation, as well as people in the so-called unresponsive wakefulness syndrome (UWS) and for people in a minimally conscious state (MCS) in the long-term rehabilitation phase F, two forms of care with a total of 70 care places are available in the "Storkower Bogen".

- Specialist inpatient care
- Outpatient care in assisted living

Our care and support concept is based on the individual wishes and personal and medical needs of our clients.

We act according to our maxim protected • required • activated.

We enable this through a highly rehabilitative approach in our care culture. This is based on the recommendations for long-term rehabilitation/ Phase F of the National Association for Rehabilitation (BAR). We also make it possible to care for and accompany our clients towards the end of their life in dignity according to recognised concepts of palliative care.

These goals can only be achieved with highly qualified employees and in an appealing, optimally equipped working and living environment.

The premises of the Centre for Ventilation and Intensive Care for inpatient care and the assisted living apartments are attractively designed. The apartments have an area of 20 m² and are comfortable, bright and friendly. All apartments and rooms are equipped with modern media and emergency call systems. The highest priority is security and privacy for our clients.

The specialization and the comprehensive range of services make the ZBI a unique facility in Berlin. We believe that high nursing and medical quality, as well as safety, can be combined with an appealing ambience, homeliness and individuality of the living environment.

Facility details

- Single apartments with more than 20 m² living space
- Wet rooms or wet area in the apartment
- Care bed, cupboard, bedside cabinet
- Therapy rooms for functional therapies
- Snoezelen rooms
- Recreation rooms on all living levels
- Bathrooms and shower baths on all living levels
- Central screen monitoring
- Screen monitoring in emergency situations
- Mobile blood gas analysis

- Compressed air and oxygen system
- Suction and oxygen units
- Nurse call system/patient call in assisted living
- AED defibrillator and emergency equipment for acute emergencies
- Central emergency power supply
- Winter garden for visitors and events

Specialist care

The outstanding feature of the concept of the ZBI Berlin is the specialisation in people in need of care with a focus on outpatient ventilation, as well as the so-called unresponsive wakefulness syndrome (UWS) and people in a minimally conscious state (MCS) Our employees are an essential guarantor for the successful implementation of our care concept. Our team consists exclusively of nurses who have passed their year 3 qualification exams from healthcare and nursing care, as well as elderly care. In addition to weekly training sessions in specialist subject areas, our employees also qualify as:

- Outpatient ventilation nurse
- Care expert for outpatient ventilation
- Care expert for people in unresponsive wakefulness syndrome/MCS
- Respiratory therapist (DGP)

In addition, we employ ergotherapy specialists, as well as respiratory and body therapists in our company. These are supported by our external therapy partners. Together, the nursing and therapy employees form a multi-professional team. Any individualised care and therapy goals are formulated and evaluated together with the clients and relatives in case discussions. Weekly therapy sessions also reflect the holistic approach to the rehabilitation process. Therapy units are effectively planned within this framework.

Based on the documented observation of 14 of our patients over several months, the good results achieved with the "Ride Java Cushion" product in terms of pressure ulcer management and ease of use, the nursing college can recommend the use of the seat cushion not only for the inpatient use, but also for the domestic environment.

6. Statement on the effectiveness of the test specimen

The test specimen was used in our facility to support pressure ulcer management for 14 patients.

The risk of pressure ulcers is assessed in our facilities based on professional nursing expertise. The Braden scale was also used for this application observation and the corresponding risk values documented.

As part of this test, the Braden scale was also used to determine the risk of pressure ulcers. In accordance with the Braden scale, all patients (residents) in our facility were at risk of developing a pressure ulcer. The following risks were documented in accordance with Braden:

- 1 patient/resident had a "medium risk"
- 13 patients had a "very high risk"

3 patients also had a pressure ulcer in the sacral region. These were pressure ulcers of categories II, III and IV (EPUAP).

In addition to using the test specimen, the patients/residents were positioned at the necessary time intervals and received individual wound care, as well as specific skin care.

Subject	Pressure ulcer localization	Wound initial situation	Wound situation after the
			three-week test period
2	Posterior region	Category III EPUAP, necrotic,	Category III EPUAP,
		redness, (depth x length x width	epithelialising, redness, (depth
		in cm) 2.4 x 6.6 x 12.4	x length x width in cm) 0.5 x
			2.5 x 4.3
3	Posterior region	Category II EPUAP, fibrinous,	Category II EPUAP,
		redness, (depth x length x width	epithelialising, redness, (depth
		in cm) 0.2 x 3.4 x 2.6	x length x width in cm) 0.0 x
			1.0 x 0.5
8	Posterior region	Category IV EPUAP,	Category IV EPUAP,
		granulating, redness, (depth x	granulating, redness, (depth x
		length x width in cm) 6.0 x 2.0	length x width in cm) 5.6 x 1.7
		x 2.0	x 1.6

The table presentation of the wound progression shows that a visible improvement in the wound condition was observed in the case of pressure ulcers.

Important components of our facility's pressure ulcer management include not only local wound treatment, but also all preventative measures such as nutrition, skin care, promoting movement, pressure distribution through individual positioning and the use of anti-decubitus mattresses and seat cushions. Consequently, the use of an anti-decubitus seat cushion is an important part of the overall nursing care concept.

From the observations and results of the test, it must be summarised that the nursing team is satisfied with the medical benefit of the test specimen, as no pressure ulcers or recurrences have occurred and the treatment of pressure ulcers has been positively supported.

7. Assessment of the benefits against the risks

At our facility, we were not able to observe any hazards arising from the product, either for patients or care staff.

The manufacturer's instructions for use are understandable and contain the most important information about product properties and materials used, as well as instructions for cleaning, maintenance and care.

8. Evaluation of desired and undesirable consequences

There were no undesirable consequences when using the product. In addition, no negative changes in the microclimate due to the occurrence of specific moisture on the skin, for example, could be observed.

9. Product description

The Ride Java seat cushion is an adjustable wheelchair cushion designed to help counteract soreness to a large extent and promote wound healing for existing pressure ulcers in the sacral region. It offers very good pressure relief for the buttock region with even pressure distribution to the buttock muscles and thighs without pressure peaks. It also enables postural correction using wedges, inserts and thigh supports. It also offers its user good sitting stability and counteracts the build-up of heat and moisture.

The cushion surrounds the pelvis with supportive contours and reduces pressure in risk

areas around the buttocks. The ventilated base cushion construction, contoured foam overlay and breathable spacer cover help to keep skin dry. The upper viscoelastic layer adapts to the body shape and helps to reduce localised pressure peaks. The contoured base helps promote and control better posture and balance.

CAM wedges can be used as accessories to further adjust the cushion. The conceptual design of the seat cushion allows easy but individual adaptation to the respective user. In combination with the breathable outer cover, a very good microclimate is created over the entire sitting period.

9.1 Changing covers

The construction of the cover allows the cover to be changed as it has a zipper.

9.2 Adaptability of the product

The product can be adjusted by the user/resident. Different wedges, medial and lateral thigh support inserts or foam inserts can be used to adapt the product to individual needs.

9.3 Usability of adjustable values by the user

The product is adjusted using wedges, thigh support inserts or foam inserts. They not only resulted in good pressure distribution, but also stabilisation in the event of posture weaknesses and balancing in the event of balance deficits. See also section 9.2.

9.4 Microclimate characteristics

There was no increased occurrence of moisture in the area of the seat surface.

9.5 Cleaning and disinfection

The cushion and cover can be easily cleaned and disinfected (see instructions for use).

9.6 Assessment of the pressure relief and distribution

The pressure relieving and redistributing effect is evaluated as positive as no new skin irritation was observed in any of the high-risk patients during the test period using the test specimen, neither in the pressure-relieved areas nor in the areas with cushion contact. The existing pressure ulcers also showed a clear healing tendency.

9.7 Ability to change positions as required

The usual seated position changes could be carried out.

9.8 Changes in the position support intervals (positioning intervals)

The position support intervals were significantly extended in almost all test subjects.

9.9 Estimation of the applied shear forces

The cushion does not adversely affect shear forces.

9.10 Changes in sitting stability

The sitting stability of the test subjects changed for the better in the majority of cases. Subjects sat more stable and slipped forward less while seated in the wheelchair using the cushion.

9.11 Changes in sitting symmetry

The sitting symmetry has also improved. Patients were much more stable and centred.

10. Overall assessment

The "Ride Java Cushion" anti-decubitus seat cushion was used successfully in our facility. Our patients/residents were satisfied with the product and the resulting sitting quality. The nursing specialists are positive towards the product after several months of testing. Handling is easy to assess for both users and carers.

The test specimen was used by high-risk patients over a period of six weeks each. Three subjects had pressure ulcers of categories II, III and IV (EPUAP) in the sacral region. Progress in wound healing was observed in all of them.

Throughout the period of use, none of our patients/residents experienced new pressurerelated skin lesions, which proved the effectiveness of the cushion.

Based on the results and experience gained during use of the cushion being tested, we assess the product as suitable for use in pressure ulcer management (pressure ulcer prevention and treatment) in accordance with the manufacturer's indications. The use of the cushion can be advocated for both domestic and inpatient use.

The content of this report or parts thereof may not be published, made accessible to the public or a group of people outside the company or used for advertising purposes without the written consent of the author.

Stade, 29/05/2019

Test plan

Type/design:

Verification of medical benefit / for products in product group 11

Project 18-123

Page 1 of 3

Art. no.: JC-2020

Product designation: 8 cm high anti-decubitus foam seat cushion

Ride Java Cushion (36x 36 cm) Art. no.: JC-1414 Ride Java Cushion (36 x 4t cm) Art. no.: JC-1416 Ride Java Cushion (38 x 38 cm) Art. no.: JC-1515 Ride Java Cushion (38 x 43 cm) Art. no.: JC-1517 Ride Java Cushion (41 x 41 cm) Art. no.: JC-1616 Art. no.: JC-1618 Ride Java Cushion (41 x 46 cm) Ride Java Cushion (41 x 51 cm) Art. no.: JC-1620 Ride Java Cushion (43 x 43 cm) Art. no.: <u>JC-1717</u> Ride Java Cushion (46 x 41 cm) Art. no.: JC-1816 Ride Java Cushion (46 x 46 cm) Art. no.: JC-1818 Ride Java Cushion (46 x 51 cm) Art. no.: JC-1820 Ride Java Cushion (51 x 41 cm) Art. no.: JC-2016 Ride Java Cushion (51 x 46 cm) Art. no.: JC-2018

incl. cover

Tested accessories: None

Manufacturer: Ride Designs a branch of Aspen Seating, LLC, 8100 SouthPark Way,

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4d,10369 Berlin

Please create a written report with the following content:

- 1. Information on product designation, type/design with item no., the tested accessories where applicable with item no., manufacturer and customer, e.g. as listed above.
- 2. Brief description of your facility/institution that describes the independence and reputation of your institution and information on the tester.
- Information on the period of use/observation Note: A minimum observation period of 6 weeks per patient is recommended.
- 4. Details of the users/subjects used for the test and the number of subjects. *Note: The test must be carried out taking into account the intended user group. We recommend including at least 10 patients in the observations.*
- 5. Information on the test environment (e.g. clinic, home setting)

 Note: The studies/tests must be carried out under general living conditions/home environment or it must be possible to apply the results to these areas.

- 6. Statement on the entire indication area specified by the manufacturer and the effectiveness of the product in combination with all tested components (e.g. covers, accessories, etc.), with:
 - 6.1. Description of the intended goal and
 - Ensuring medical treatment (decubitus therapy) and/or
 - Prevention of a disability (decubitus prophylaxis) and/or
 - Treatment of a disability (e.g. improving freedom of movement by changing position)
 - 6.2. Clear description of the goal(s) achieved and outlined under 6.1., along with reason why the tester believes that the product is suitable for use in the home environment with the indications made by the manufacturer, complete with clinical endpoints, the pressure ulcer incidences (occurrence of pressure ulcer stage II or higher) to verify the effectiveness at preventing pressure ulcers, or, for products intended for pressure ulcer treatment, the healing rate (specifying the measurement method use) of an occurring pressure ulcer. Here, the indication area used must be clearly formulated and outlined based on standardised risk scales the Braden scale in the case of geriatric patients and other approved scales for other patient groups; if the product is intended for pressure ulcer treatment, a description of the relevant pressure ulcer stages must also be provided in accordance with EPUAP.
- 7. An assessment of the benefits against the risks.
- 8. An assessment of the desired and undesirable consequences ("outcomes")
- 9. Description of the product with statement on the following points:
 - 9.1. Can the cover be changed by the user?
 - 9.2. Can the product be adapted to the individual weight of the patient and the relevant load situation (manually or automatically)? Or can the product be selected based on the weight of the patient?
 - 9.3. Can all individually adjustable values be operated by the user?
 - 9.4. Observable microclimatic properties
 - 9.5. Further information on the handling of the product where applicable, e.g. cleaning/disinfection, comfort for the user/patient and care staff.
 - 9.6. Assessment of the pressure relief and distribution
 - 9.7. Ability to change positions as required
 - 9.8. Changing the position support intervals (positioning intervals)
 - 9.9. Estimation of the applied shear forces
 - 9.10. Changes in sitting stability
 - 9.11. Changes in sitting symmetry

10. Overall assessment

The report can be written in German or English and pictures and evaluation tables may be added where applicable.

Confirmation of the test facility /expert

I hereby confirm that the assessment of the medical benefit for the product listed above was created according to the test plan provided.

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Das Leben pflegen...

0 Sandersweg 15 in 21680 Stade 2 04141- 535 9810 ^ 04141 -5359830 S www.johannisheim-stade.de

Johannisheim e.V. - Sandersweg 15 in 21680 Stade

BEO MedConsulting Berlin GmbH attn: Mr. Friedrich Helmholtzstr. 2

Nursing services Sandra Idem S 04141-5359828 (E3 sandra.idem@iohannisheim-stade.de

22/03/2019

10587 Berlin - Charlottenburg

Second supplementary evidence of the medical benefit of the Ride Forward Cushion antidecubitus foam seat cushion to assist in the prevention and treatment of pressure ulcers by Ride Designs, USA

Dear Mr. Friedrich.

I am pleased to send you today my second and supplementary statement on the proof of the medical benefit of the "Ride Forward Cushion" anti-decubitus seat cushion. As requested, we have used the cushion on another four residents in our facility.

The second test phase confirms the results of the previous assessment.

Should you have any questions about the results presented, please contact me.

I hope I can help you with my explanations.

Best regards,

Sandra Idem

- Nurse & nursing services -

1. The "Ride Forward Cushion" anti-decubitus seat cushion incl. protective cover was evaluated again. The product is intended to assist in the prevention and treatment of pressure ulcers in accordance with the manufacturer's product documentation. It is offered in the following dimensions:

Ride Forward Cushion (25 x 25 cm) Art. no.: <u>FCD-1010</u> Ride Forward Cushion (25 x 30 cm) Art. no.: <u>FCD-1012</u> Ride Forward Cushion (30 x 30 cm) Art. no.: <u>FCD-1212</u> Ride Forward Cushion (30 x 36 cm) Art. no.: <u>FCD-1214</u> Ride Forward Cushion (30 x 41 cm) Art. no.: FCD-1216 Ride Forward Cushion (36 x 36 cm) Art. no.: FCD-1414 Ride Forward Cushion (36 x 41 cm) Art. no.: FCD-1416 Ride Forward Cushion (36 x 46 cm) Art. no.: FCD-1418 Ride Forward Cushion (38 x 38 cm) Art. no.: FCD-1515 Ride Forward Cushion (38 x 43 cm! Art. no.: FCD-1517 Ride Forward Cushion (41 x 41 cm) Art. no.: FCD-1616 Ride Forward Cushion (41 x 46 cm) Art. no.: <u>FCD-1618</u> Ride Forward Cushion (41 x 51 cm) Art. no.: <u>FCD-1620</u> Ride Forward Cushion (43 x 43 cm) Art. no.: FCD-1717 Ride Forward Cushion (46 x 41 cm) Art. no.: FCD-1816 Ride Forward Cushion (46 x 46 cm) Art. no.: <u>FCD-1818</u> Ride Forward Cushion (46 x 51 cm) Art. no.: FCD-1820 Ride Forward Cushion (51 x 41 cm) Art. no.: <u>FCD-2016</u> Ride Forward Cushion (51 x 46 cm) Art. no.: <u>FCD-2018</u> Ride Forward Cushion (51 x 51 cm) Art. no.: <u>FCD-2020</u>

incl. inner and outer

protective cover

The scope of delivery includes the cushion with an inner and outer protective cover.

2. Information on the tester and institution

The Ride Forward Cushion anti-decubitus seat cushion was tested at the Johannisheim care home in Stade. The facility provides individual inpatient care for the elderly and the disabled (at all levels of care) in various living areas. The individual wishes and needs of guests and residents are at the heart of the nursing care.

Ms. Sandra Idem acts as the expert. She is a registered nurse and provider of nursing services. Ms Idem has been testing anti-decubitus systems for many years.

3. Information on the period of use/observation

The anti-decubitus seat cushion was tested for the first time from 15 August 2018 to 8 November 2018 in our facility with six residents over a period of six weeks each.

Furthermore, the cushion was used by four test persons over a six-week period from 21 January 2019 to 4 March 2019. The corresponding results of this second phase of the test are presented below.

4. Description of test subjects and application history

Subjects were between 32 and 86 years of age and had different underlying medical conditions. These were predominantly typical, predominantly age-related neurological, internal and neurodegenerative diseases. The risk of pressure ulcers was assessed according to the national expert standard for pressure ulcer prevention based on nursing expertise. In addition, the Braden scale was used as a standardised assessment.

Residents	07	08
Test location		Stade
Test period	21 January 20	19 to 4 March 2019
Gender	Male	Female
Age	86	32
Height in cm	175	172
Weight in kg	71	69
Diagnoses, special features of the care situation	Dementia	Premature birth (week of pregnancy 30), brain damage, spastic paresis arms, immobile
Level of care	• III (request for upgrade in progress)	• IV
Pressure ulcer risk assessment through nursing expertise	Medium risk	High risk
Braden-scale assessment of pressure ulcer risk	General risk	High risk

Special risk factors	 Resident is advanced in years The resident is suffering from malnutrition Resident is lying on a high-level pressure ulcer There are no movement restrictions Pain symptoms Dementia with aggressive characteristics No sensory stimuli perceived in the sitting area. The resident also does not consciously relieve the load on the sitting area 	 There are severe movement restrictions Paralysis, contractures and spasticity are present Sensory perception is limited No perception of sensory stimuli in the sitting area Consequently, no conscious relief of the sitting area There is no stability while sitting Sense of balance is limited The resident cannot change sitting position adequately without help The resident sits asymmetrically
Length of use of test specimen per day	 The cushion was used for an average of 8 hours during the day. The longest use intervals of the cushion without interruption were 4 hours. 	The cushion is used continuously for around 10 hours a day.
Pressure ulcer present	Two pressure ulcers in the sacral region	Ischium, left
	III and IV	I
Healing progress	Clearly visible	Clearly visible
Has there been a pressure ulcer in the past?	Tot clear from the documents available to as.	No
Estimation of pressure relief with the aid		The aid has helped in our opinion, as the wound improved.
Positioning	· · · · · · · · · · · · · · · · · · ·	The resident was positioned at the required time intervals

Residents	09	10
Test location	Stade	
Test period	21 January 2019 to 4 March 2019	
Gender	Male	Female
Age	83	65
Height in cm	182	168
Weight in kg	79.5	72
Diagnoses, special features of the care situation	Post apoplexy with right hemiparesis	MS Post apoplexy Post cancer
Level of care	• IV	• IV
Pressure ulcer risk assessment through nursing	High risk	Medium risk

expertise		
Braden-scale	Medium risk	Low risk
	iviedium risk	LOW IISK
assessment of		
pressure ulcer risk		
Special risk factors	 Resident is advanced in years Bed confinement Significant movement restrictions Paralysis Spasticity Limited sensory faculties No sensory stimuli perceived in the sitting area No conscious relief of the sitting area No stability while sitting Limited balance Inadequate change of sitting position 	Moderate movement restrictions Paralysis Limited sensory faculties Partial conscious relief of the sitting area
Length of use of test specimen per day	 Continuously for 2 hours at a time, 2 hours daily in total 	 Continuously for 4 hours at a time, approx. 8 hours daily in total
Pressure ulcer	Ischium, right	Sacral region, left
present	, ,	
	II	Ī
Healing progress	Significant improvement	Positive - the ulcer has healed
Has there been a	No	No
pressure ulcer in the past?		
Estimation of	Pressure redistribution can be assessed as	Good - the ulcer healed within the first three weeks
pressure relief with	positive, as the wound showed a clear healing	
the aid	tendency	
Positioning	The resident can be positioned in the usual way	The resident can be positioned in the usual way

5. Information on the test equipment

The "Ride Forward Cushion" anti-decubitus seat cushion was tested at our facility in Stade.

Our facility

Our care team at Johannisheim nurtures and cares for people with mental and physical complaints at all levels of care. Johannisheim has capacity for 124 residents. 90 single and 17 double rooms of different room size(s) are available. Our services cover short and long term care and inclusive dementia care. Johannisheim also cares for residents who have suffered from a stroke, Parkinson's disease, apallic conditions, Alzheimer's disease and MS, etc.

Our philosophy

Johannisheim is a social welfare institution based on the Christian ideology. This affects the basic social diaconal attitude and actions in all areas of the facility's life and work.

The job is a vocation for us

Everyone has an unrestricted right to a dignified life - even at an advanced age. This conviction determines how we deliver our service to people in care and support. We want to contribute to humanity in our society through expert support and personal care - by helping and listening.

Johannisheim sees itself as a "company of humanity"

- Our work is characterized by appreciation for people in need of care.
- We want to combine our professional competence with personal care.
- We see our cooperation as an important asset for our residents.
- The satisfaction of the people entrusted to us and who we work to help is an important benchmark for us.
- Good cooperation is an essential foundation for our welfare services.

Feedback is a development aid for us

The senior advisory board freely chosen in our company is also our in-house supervisory body and monitors all services and activities in the business, such as care, catering, supervision, accommodation, care contracts, construction measures, investments, projects, etc., but also the socially responsible treatment of the residents entrusted to us. Our senior advisory board, as well as relatives, guests and visitors, improve our quality of work and enrich our development through critical feedback of a positive and negative nature.

6. Statement on the effectiveness of the test specimen

The test specimen was used in our facility to support the management of pressure ulcers for ten residents. This information relates to the first report with six residents from December 2018 and the second report with a further four test subjects from today.

The risk of pressure ulcers is assessed in our facilities based on professional nursing expertise. In line with this assessment of the risk of pressure ulcers, six residents posed a "high risk" and four a "medium risk".

As part of this test, the Braden scale was also used to determine the risk of pressure ulcers. In accordance with the Braden scale, all residents in our facility were at risk of developing a pressure ulcer. The following risks were documented in accordance with Braden:

- 2 resident had a "general risk"
- 2 resident had a "low risk"
- 4 residents had a "medium risk"
- 2 resident had a "very high risk"

In addition, nine residents had one pressure ulcer in the sacral region and ischiums - five pressure ulcers in category I (EPUAP), three in category II (EPUAP), another in category III (EPUAP) and one in category IV (EPUAP).

In addition to using the test specimen, the residents were positioned at the necessary time intervals and received individual wound care, as well as specific skin care.

Subject	Pressure ulcer localization	Wound initial situation	Wound situation after the three-week test period
1	Posterior region	Category II EPUAP, fibrinous, no redness (depth x length x width in cm) 0.2 x 2.0 x 1.5	Category II EPUAP, granulating, no redness (depth x length x width in cm) 0.1 x 1.5 x 1.3
2	Posterior region	Category I EPUAP, (depth x length x width in cm) 0.0 x 1.0 x 2.0	Category I EPUAP, (depth x length x width in cm) 0.0 x 1.0 x 1.0
3	Posterior region	Category II EPUAP, fibrinous, signs of wound infection (depth x length x width in cm) 0.1 x 2.0 x 1.0	fibrinous, signs of wound
4	Posterior region	Category I EPUAP, (depth x length x width in cm) 0.0 x 2.0 x 2.0	Healed
5	Posterior region	Category I EPUAP, (depth x length x width in cm) 0.0 x 2.0 x 2.0	Healed
7	Sacral region (A)	Category III EPUAP, fibrinous, redness, (depth x length x width in cm) 0.3 x 2.0 x 1.5	
7	Sacral region (B)	Category III EPUAP, fibrinous, redness, (depth x length x width in cm) 0.4 x 2.0 x 2.0	
7	Trochanter right	Category IV EPUAP, necrotic, redness, (depth x length x width in cm) 0.5 x 9.0 x 8.0	

8	Ischium, left	length x width in cm) 0.0 x 1.0	Category I EPUAP, (depth x length x width in cm) 0.0 x 0.5 x 0.5
9	Ischium, right	,	
10	Sacral region, left	length x width in cm) 0.0 x 2.0	The pressure ulcer healed within the third week of testing.

The table presentation of the wound progression shows that a visible improvement in the wound condition was observed in the case of pressure ulcers or three of the pressure ulcers healed.

Important components of our facility's pressure ulcer management include not only local wound treatment, but also all preventative measures such as nutrition, skin care, promoting movement, pressure distribution through individual positioning and the use of anti-decubitus mattresses and seat cushions. Consequently, the use of an anti-decubitus seat cushion is an important part of the overall nursing care concept.

From the observations and results of the test, it must be summarised that the nursing team is satisfied with the medical benefit of the test specimen, as no pressure ulcers or recurrences have occurred and the treatment of pressure ulcers has been positively supported.

7. Assessment of the benefits against the risks

At our facility, we were not able to observe any hazards arising from the product, either for patients or care staff.

The manufacturer's instructions for use are easily understandable and contain the most important information about product properties and materials used, as well as instructions for cleaning, maintenance and care.

8. Evaluation of desired and undesirable consequences

There were no undesirable consequences when using the product. Furthermore, for example, no negative changes in the microclimate were observed due to the occurrence of specific moisture on the skin, although, during the first test (in August and September), there were sometimes high summer outside temperatures.

These observations were confirmed in the second test phase. Occurrence of above-average skin moisture was also not observed during this test period.

9. Product description

The "Ride Forward Cushion" anti-decubitus wheelchair seat cushion is a medical device intended for the prevention, treatment or alleviation of pressure ulcers in the sacral/ischium region.

The cushion consists of:

- The actual seat cushion core
- An inner cover
- An outer cover.

The seat cushion core consists of two layers. The upper part is made of viscoelastic foam, the lower part of polyurethane foam. The foam core is anatomically preformed.

The seat cushion core is covered by a first cover made of polyether. It is equipped with a zipper and protects the core from the penetration of liquids.

The outer cover, which is also removable, is made of air-permeable fabrics that are intended to counteract the increased occurrence of moisture on the skin. The underside is also slip-resistant, preventing the cushion from "slipping" on the support (e.g. the seat of a wheelchair).

9.1 Changing covers

The construction of the two covers allows them to be changed as they both have a zipper.

9.2 Adaptability of the product

The product cannot be adjusted by the user.

To use the cushion, simply place it on the seat base. The "Ride Forward Cushion" anti-decubitus seat cushion has a maximum user weight. It can be used for patients weighing 5 to 113 kg.

9.3 Usability of adjustable values by the user The product cannot be changed by the user.

9.4 Microclimate characteristics

There was no increased occurrence of moisture in the area of the seat surface in both test phases.

9.5 Cleaning and disinfection

It is possible to clean and disinfect the cushion as well as the covers (see instructions for use).

9.6 Assessment of the pressure relief

The pressure-distributing effect is evaluated as positive, as no new skin irritation was observed in any of the at-risk patients during the test period for the test specimen. The existing pressure ulcers also showed a healing tendency or healed.

9.7 Ability to change positions as required

The usual seated position changes could be carried out.

9.8 Changes in the position support intervals (positioning intervals)

The position support intervals remained practically unchanged for all subjects during the test phase.

9.9 Estimation of the applied shear forces

The cushion does not adversely affect shear forces.

9.10 Changes in sitting stability

The sitting stability of the test subjects has not changed negatively.

9.11 Changes in sitting symmetry

The sitting symmetry of the test subjects has not changed negatively.

10. Overall assessment

The "Ride Forward Cushion" anti-decubitus seat cushion was successfully tested in our care facilities. Our residents were satisfied with the product and the resulting sitting quality. The nursing specialists are positive towards the product after several months of testing. Handling is easy to assess for both users and carers.

The test specimen was used by high-risk patients over a period of six weeks each. Nine subjects had pressure ulcers of categories I, II, III and IV (EPUAP) in the sacral region or in the ischium area. Progress in wound healing was observed in all of them. Three of the pressure ulcers healed.

Throughout the period of use, none of our residents experienced new pressure-related skin lesions, which proved the effectiveness of the cushion.

Based on the results and experience gained during use of the cushion being tested, we assess the product as suitable for use in pressure ulcer management (pressure ulcer prevention and treatment) in accordance with the manufacturer's indications. The use of the cushion can be advocated for both domestic and inpatient use.

The content of this report or parts thereof may not be published, made accessible to the public or a group of people outside the company or used for advertising purposes without the written consent of the author.

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